Intended use

Easicult Combi slides are intended for monitoring of microbial contamination in various industrial environments. The test can be performed on-site, or the slides can be used as convenient transport media for samples.

The slide is covered on one side with TTC Agar (yellowish), which supports growth of most common bacteria, and on the other side with Rose-Bengal Agar (pink), which supports growth of fungi.

The main significance of the test is that elevation of total microbial counts can be detected.

Contents of the kit

Easicult Combi	Cat. Nos. 67987, 05984
Test slides	10 pcs
Labels	10 pcs
Instructions for use	1 pc

Typical formulation

TTC Agar	Rose-Bengal Agar
Tryptone	Peptone
Soy peptone	Dextrose
Disodium succinate	Potassium hydrogen phosphate
TTC solution	Magnesium sulphate
Agar-agar	Sodium chloride
Water	Rose Bengal
	Sodium hydroxide
	Gentamycin sulphate
	Chloramphenicol
	Agar-agar
	Water

Warnings and precautions

Do not use the product beyond the expiry date marked on the kit.

Do not touch the unused growth media.
Store in upright position.

Do not use the slide if you notice

- discoloration or dehydration of the growth medium
- detachment of the growth media from the plastic slide
- evidence of microbial growth.

Because any growth on the Easicult Combi slide may be pathogenic, do not touch the used slide.

Storage

Store Easicult Combi at room temperature (18...25°C / 64...77°F) protected from draught, temperature fluctuations and light sources. Avoid storage near heat-generating appliances. Do not allow to freeze. The expiry date (year-month-date) is marked on the box and on the cap of each slide.

Sampling and procedure (Fig. 1–5)

To avoid contamination, the growth medium should not come into contact with any other material than the one to be tested. On the other hand, it is important that the growth medium makes full contact with the material to be tested.

Viscous fluids and fluids with high bacterial content (>10⁷ CFU/ml)

If the viscosity or bacterial content of the sample is high, the sample should be diluted. For dilution, put 100 or 1000 ml of drinkable tap water in a clean, well-rinsed and dried bottle with a cap. The bacterial content of the water for dilution should not exceed 100 CFU/ml. Before filling the bottle, let the water run for 5 minutes or boil it for 15 minutes and then let it cool. Using a clean (disposable) pipette, add 1 ml of sample to the bottle. Close and mix thoroughly by shaking the bottle about 30 times. Dip the slide in the dilution and proceed as described for fluid samples.

Fluid samples

- 1 Unscrew the tube and withdraw the slide without touching the agar surfaces.
- 2 Dip the slide in the liquid. Alternatively, wet the slide under a running stream of the liquid or spray the liquid on the slide. If the liquid is under pressure, the slide must be handled carefully to avoid unfastening of the agar. If the sample is in a container, mix the contents and dip the slide in it. Both agar sides should get completely wet. The slide must be in contact with the fluid for 5 to 10 seconds.
- 3 Allow excess fluid to drain off the slide. Blot the last drops from the lower end of the slide on absorbent paper.

- 4 After sampling, screw the slide tightly back into the tube. Fill in the label and affix it to the tube.

- 5 Incubate the slide at 27...30°C (80...86°F) for 24–48 hours for detection of bacteria. Yeasts and moulds need three days' incubation before the result can be read. If the incubation is carried out at room temperature, the incubation times are 2–4 days and 4–7 days, respectively.

If the normal temperature for the fluid tested differs substantially from the above incubation temperatures, it may result in slow growth of bacteria. An incubation temperature closer to the normal temperature of the fluid is then recommended.

Interpretation of results (Fig. 6)

Cautiously remove the slide from its tube after incubation and determine the microbial count (number of colony forming units, CFU) by comparing the density of growth on the slide with the model chart. Note the separate model charts for TTC and Rose-Bengal.

If the sample was diluted, the dilution factor must be taken into account in the evaluation. For example, if a dilution of 1 + 100 ml (1 ml of sample in 100 ml of water) shows a density of 10⁶ CFU/ml, the actual result for the sample is 10⁸ CFU/ml.

Rose-Bengal medium, determination of fungi

Growth on Rose-Bengal Agar may consist exclusively of either yeasts or moulds, or mixed growth of both. Mould colonies are soft and fluffy and usually pale, green or black in colour. Yeasts usually grow in dome-shaped colonies but may sometimes be flat and dry. Since mould colonies may originate from fragments of mycelium or from individual spores, the model chart is not quantitative. The chart indicates whether the contamination is slight (+), moderate (++) or heavy (+++). Colonies can be removed from the slides and examined under microscope if needed. Fungal contamination is sometimes apparent to the naked eye as a coating on the fluid surface.

TTC Agar, determination of total bacterial count

Most aerobic bacteria grow on TTC Agar as red colonies. Moulds and yeasts may also grow slowly on this medium. Even though the bacterial growth is almost always in the form of red colonies, any colourless colonies should also be taken into account when the density is estimated. In cases where large colonies are present, it should be borne in mind that colony density, not the size of individual colonies, should be considered.

If the bacterial count is very high (more than 10⁷ CFU/ml), the growth is confluent. This can appear as a uniformly red surface layer. Very rarely there is totally colourless growth. It is advisable to compare slides exhibiting an apparently uniform surface with an unused slide to avoid misinterpretation.

As no universally applicable limits exist, limit values have to be determined by experience. For bacterial contamination in coolants, the following guide may be used:

Bacteria		
CFU/ml	Contamination	
< 10 ⁴	slight	usually no problems ¹
10 ⁴ –10 ⁶	moderate	
> 10 ⁶	heavy	not acceptable ¹

Limitations of the method

If the bacterial count exceeds 10⁷ CFU/ml, or the viscosity is high, the sample should be diluted.

Very rarely the bacteria grow on the TTC medium as colourless colonies.

The reliable lower detection limit for bacteria is 10³ CFU/ml. The growth of some coccoid bacteria may be weakened by TTC.

Disposal

- Dispose of contents according to national and local law.
- All used components should be handled and disposed of as potentially pathogenic material.
- Materials of the components:
Paper: Instructions for use, patient labels
Cardboard: Kit box
Plastic: Tubes, caps and dipslides
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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Easicult® is a registered trademark of Aidian Oy.